

## Product datasheet for **TP723239**

### **IL5 (NM\_021834) Rat Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Purified recombinant protein of Rat interleukin 5 (IL5).
<b>Species:</b>	Rat
<b>Expression Host:</b>	E. coli
<b>Expression cDNA Clone or AA Sequence:</b>	MEIPMSTVVK ETLIQLSTHR ALLTSNETMR LPVPTHKNHQ LCIGEIFQGL DILKNQTVRG GTVEILFQNL SLIKKYIDGQ KEKCGEERRK TRHFLDYLQE FLGVMSTEWA MEV
<b>Tag:</b>	Tag Free
<b>Predicted MW:</b>	26 kDa
<b>Concentration:</b>	lot specific
<b>Purity:</b>	>95% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	Lyophilized from a 0.2 µM filtered solution of 20mM phosphate buffer, 100mM NaCl, pH 7.2
<b>Bioactivity:</b>	ED50 was determined by the dose-dependent stimulation of the proliferation of human TF-1 cells is > 0.5 ng/ml, corresponding to a specific activity of > 2 x 10 <sup>6</sup> units/mg.
<b>Endotoxin:</b>	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
<b>RefSeq:</b>	<a href="#">NP_068606</a>
<b>Locus ID:</b>	24497
<b>UniProt ID:</b>	<a href="#">Q9R2C9</a>
<b>RefSeq Size:</b>	399
<b>Cytogenetics:</b>	10q22
<b>RefSeq ORF:</b>	396
<b>Summary:</b>	cytokine with B-cell growth factor activity [RGD, Feb 2006]



[View online »](#)

**Product images:**

