

Name: Purified recombinant protein of Rat interleukin 6 (Il6).

Catalog: TP723242

Product Data Sheet - PROTEIN

Gene Name: interleukin 6

GeneBank accession: NM_012589

Description: Purified recombinant protein of Rat interleukin 6 (Il6). (TP723242)

Amount: 10 ug

Buffer and Storage: Lyophilized from a 0.2 µM filtered solution of 20mM phosphate buffer, 100mM NaCl, pH 7.2

Source: Recombinant protein was produced with
MFPTSQVRRGDFTEDTTHNRPVYTTSQVGGILITYVLRILEMRKELCNGNSDCMNSDDA
LSENNLKLPEIQRNDGCFQTGYNQEICLLKICSGLLEFRFYLFVKNLQDNKKDKARVIQS
NTETLVHIFKQEIKDSYKIVLPTPTSNALLMEKLESQKEWLRTKTIQLILKALEEFLKVTMRST
RQT

Protein Accn: NP_036721

Gene Synonym: IL6|Ifnb2

Predicted molecular weight: 21.7 kDa

Tags: Tag-free

Purity: >95% as determined by SDS-PAGE and Coomassie blue staining

Concentration: Resuspend the protein in the desired concentration in proper buffer

* More data may be available on our website: [TP723242](https://www.origene.com/TP723242)