

Name: Purified recombinant protein of Mouse chemokine (C-C motif) ligand 7 (Ccl7). Product Data Sheet - PROTEIN	Catalog: TP723285
--	--------------------------

Gene Name: chemokine (C-C motif) ligand 7

GeneBank accession: NM_013654

Description: Purified recombinant protein of Mouse chemokine (C-C motif) ligand 7 (Ccl7). (TP723285)

Amount: 10 ug

Buffer and Storage: Lyophilized from a 0.2 µM filtered solution of 20mM phosphate buffer, 100mM NaCl, pH 7.2

Source: Recombinant protein was produced with
QPDGPNASTCCYVKKQKIPKRNLSYRRITSSRCPWEAVIFKTKKGMEVCAEAHQKWVEE
AIAYLDMKTPTPKP

Protein Accn: NP_038682

Gene Synonym: MCP-3|Scya7|fic|marc|mcp3

Predicted molecular weight: 8.5 kDa

Tags: Tag-free

Purity: >95% as determined by SDS-PAGE and Coomassie blue staining

Concentration: Resuspend the protein in the desired concentration in proper buffer

* More data may be available on our website: [TP723285](#)