

# **Product datasheet for TP723748**

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## II10 (NM\_010548) Mouse Recombinant Protein

#### **Product data:**

**Product Type:** Recombinant Proteins

**Description:** Purified recombinant protein of Mouse interleukin 10 (II10)

Species: Mouse Expression Host: E. coli

**Expression cDNA Clone** Mouse IL-10, the region of Ser19-Ser178, from gene Accession# NM\_010548, with a Cys149Tyr

or AA Sequence: substitution

Tag: Tag Free
Predicted MW: 18.9 kDa
Concentration: lot specific

**Purity:** Purity is >98%, as determined by Coomassie stained SDS-PAGE.

Buffer: 20 mM Tris pH 7.5

**Bioactivity:** The ED50 is 1- 2 ng/ml, corresponding to a specific activity of 1- 0.5 x 10<sup>6</sup> units/mg,

determined by the dose dependent stimulation MC/9 cell proliferation.

**Endotoxin:** Less than 0.01 ng per ug protein as determined by the LAL method

Storage: Store at -80°C.

Stability: Unopened vial can be stored between 2°C and 8°C for up to 2 weeks, at -20°C for up to 6

months, or at -70°C or below until the expiration date. Aliquots can be stored between 2°C and 8°C for up to one week and stored at -20°C or colder for up to 3 months. Avoid repeated

freeze/thaw cycles.

**RefSeq:** NP 034678

**Locus ID:** 16153

UniProt ID: P18893, Q3U879

RefSeq Size: 1306

Cytogenetics: 1 56.89 cM

RefSeq ORF: 534

Synonyms: CSIF; If2a; IL-; Il-10





### **Summary:**

This gene encodes an anti-inflammatory cytokine that is a member of the class-2 cytokine family. The encoded protein is secreted by cells of both the innate and adaptive immune systems and is crucial for limiting the immune response to a broad range of pathogens. It also has been shown to suppress autoimmune responses. This protein mediates it's immunosuppressive signal through a specific interleukin 10 receptor complex. Aberrant functioning of this gene is associated with numerous immune disorders including graft-versus-host disease, and increased susceptibility to HIV-1 infection and rheumatoid arthritis. [provided by RefSeq, Sep 2015]

# **Product images:**

