

Product datasheet for **TP723913**

SARS-CoV-2 Spike ECD Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	SARS-CoV-2 Spike ECD Recombinant Protein composed of 1209 amino acids from Spike Protein S1, S2 and ECD domains, with His-tag
Species:	SARS-CoV-2
Expression Host:	Insect
Tag:	His-tag
Predicted MW:	134 kDa (His-tag)
Purity:	> 90%
Endotoxin:	10 mM Phosphate Buffered Saline, pH 7.4
Storage:	Store at 2-8°C for short Term (<2 weeks). Long Term: -20°C. For long term storage, aliquot, and store at -20°C.
Stability:	1 year
Locus ID:	43740568



[View online »](#)