

Product datasheet for TP723921


SARS-CoV-2 Spike S2 ECD Recombinant Protein with His tag

Product data:

Product Type:	Recombinant Proteins
Description:	SARS-CoV-2 Spike S2 ECD Recombinant Protein with His tag
Species:	SARS-CoV-2
Expression Host:	Insect
Expression cDNA Clone or AA Sequence:	SVASQSIAY TMSLGAENSV AYSNNSIAIP TNFTISVTTE ILPVSMTKTS VDCTMYICGD STECSNLLLQ YGSFCTQLNR ALTGIAVEQD KNTQEVEFAQV KQIYKTPPIK DFGGFNFSQI LPDPSKPSKR SFIEDLLFNK VTLADAGFIK QYGDCLGDIA ARDLICAQKF NGLTVLPPLL TDEMIAQYTS ALLAGTITSG WTFGAGAALQ IPFAMQMAYR FNGIGVTQNV LYENQKLIAN QFNSAIGKIQ DSLSSTASAL GKLQDVVNQN AQALNTLVKQ LSSNFGAISS VLNDILSRLD KVEAEVQIDR LITGRLQSLQ TYVTQQLIRA AEIRASANLA ATKMSECVLG QSKRVDFCGK GYHLMSFPQS APHGCVVFLHV TYVPAQEKNF TTAPAICHDG KAHFPREGVF VSNGTHWFVT QRNFYEPQII TTDNTFVSGN CDVVIGIVNN TVYDPLQPEL DSFKEELDKY FKNHTSPDVD LGDISGINAS VVNIQKEIDR LNEVAKNLNE SLIDLQELGK YEQYIKWP
Tag:	C-His
Purity:	> 90% by SDS-PAGE.
Buffer:	Lyophilized from a 0.22 µm filtered solution of sterile 20 mM PB, 300 mM NaCl, pH 7.0.
Endotoxin:	< 1.0 EU/µg of the protein by LAL method.
Protein Description:	Recombinant 2019-nCoV Spike S2 ECD Protein is produced by Baculovirus-Insect Cells expression system. The target protein is expressed with sequence (Ser686-Pro1213) of 2019-ncov Spike S2 ECD (Accession #YP_009724390.1) fused with a His tag at the C-terminus.
Storage:	Store at -20°C to -80°C.
Stability:	3 months
Locus ID:	43740568



[View online »](#)

Product images:

Recombinant SARS-COV-2 Spike S2 ECD Protein with His tag was determined by SDS-PAGE with Coomassie Blue, showing a band at 70 kDa.