

Product datasheet for **TP724234**

Mouse ADAM9 Protein, His Tag

Product data:

Product Type:	Recombinant Proteins
Description:	Mouse ADAM9 Protein, His Tag
Expression Host:	HEK293
Tag:	C-6×His
Predicted MW:	The protein has a predicted molecular mass of 74.3 kDa after removal of the signal peptide. The apparent molecular mass of mADAM9-His is approximately 55-70 kDa due to glycosylation.
Purity:	The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.
Reconstitution Method:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Stability:	12 months from date of despatch
Summary:	Cleaves and releases a number of molecules with important roles in tumorigenesis and angiogenesis, such as TEK, KDR, EPHB4, CD40, VCAM1 and CDH5 (PubMed:19273593). May mediate cell-cell, cell-matrix interactions and regulate the motility of cells via interactions with integrins (PubMed:10825303).[UniProtKB/Swiss-Prot Function]



[View online »](#)