

Product datasheet for **TP724249**

Mouse TENM4 Protein, His Tag

Product data:

Product Type:	Recombinant Proteins
Description:	Mouse TENM4 Protein, His Tag
Expression Host:	HEK293
Tag:	N-6×His
Predicted MW:	The protein has a predicted molecular mass of 269.0 kDa after removal of the signal peptide.
Purity:	The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.
Reconstitution Method:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Stability:	12 months from date of despatch
Summary:	Involved in neural development, regulating the establishment of proper connectivity within the nervous system. Plays a role in the establishment of the anterior-posterior axis during gastrulation. Regulates the differentiation and cellular process formation of oligodendrocytes and myelination of small-diameter axons in the central nervous system (CNS). Promotes activation of focal adhesion kinase. May function as a cellular signal transducer. [UniProtKB/Swiss-Prot Function]



[View online »](#)