

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TP724269

Human HPSE Protein, His Tag

Product data:

| Product Type: | Recombinant Proteins |
|------------------------|---|
| Description: | Human HPSE Protein, His Tag |
| Expression Host: | HEK293 |
| Tag: | C-6×His |
| Predicted MW: | The protein has a predicted molecular mass of 58.4 kDa after removal of the signal peptide. The apparent molecular mass of HPSE-His is approximately 55-70 kDa due to glycosylation. |
| Purity: | The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining. |
| Reconstitution Method: | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. |
| Storage: | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Stability: | 12 months from date of despatch |
| Summary: | Heparan sulfate proteoglycans are major components of the basement membrane and extracellular matrix. The protein encoded by this gene is an enzyme that cleaves heparan sulfate proteoglycans to permit cell movement through remodeling of the extracellular matrix. In addition, this cleavage can release bioactive molecules from the extracellular matrix. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011] |



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US