

## Product datasheet for **TP724302**

### Human LIV-1 Protein, His Tag

#### Product data:

|                               |   |
|-------------------------------|---|
| <b>Product Type:</b>          | Recombinant Proteins  |
| <b>Description:</b>           | Human LIV-1 Protein, His Tag  |
| <b>Expression Host:</b>       | HEK293  |
| <b>Tag:</b>                   | C-6×His   |
| <b>Predicted MW:</b>          | The protein has a predicted molecular mass of 34.4 kDa after removal of the signal peptide. The apparent molecular mass of LIV-1-His is approximately 55-70 kDa due to glycosylation.   |
| <b>Purity:</b>                | The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.  |
| <b>Reconstitution Method:</b> | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.  |
| <b>Storage:</b>               | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.   |
| <b>Stability:</b>             | 12 months from date of despatch   |
| <b>Summary:</b>               | Zinc is an essential cofactor for hundreds of enzymes. It is involved in protein, nucleic acid, carbohydrate, and lipid metabolism, as well as in the control of gene transcription, growth, development, and differentiation. SLC39A6 belongs to a subfamily of proteins that show structural characteristics of zinc transporters (Taylor and Nicholson, 2003 [PubMed 12659941]).[supplied by OMIM, Mar 2008] |



[View online »](#)