

Product datasheet for **TP724349**

Human GPR87 Protein, mFc Tag

Product data:

Product Type:	Recombinant Proteins
Description:	Human GPR87 Protein, mFc Tag
Expression Host:	HEK293
Tag:	C-Mouse Fc
Predicted MW:	The protein has a predicted molecular mass of 31.3 kDa after removal of the signal peptide. The apparent molecular mass of GPR87-mFc is approximately 35-55 kDa due to glycosylation.
Purity:	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Reconstitution Method:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Stability:	12 months from date of despatch
Summary:	This gene encodes a G protein-coupled receptor and is located in a cluster of G protein-couple receptor genes on chromosome 3. The encoded protein has been shown to be overexpressed in lung squamous cell carcinoma (PMID:18057535) and regulated by p53 (PMID:19602589). [provided by RefSeq, Nov 2011]



[View online »](#)