

Product datasheet for **TP724434**

Human CLDN3(144-156) Protein, mFc Tag

Product data:

Product Type:	Recombinant Proteins
Description:	Human CLDN3(144-156) Protein, mFc Tag
Expression Host:	HEK293
Tag:	C-mouse Fc
Predicted MW:	The protein has a predicted molecular mass of 27.8 kDa after removal of the signal peptide. The apparent molecular mass of CLDN3(144-156)-mFc is approximately 25-35 kDa due to glycosylation.
Purity:	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Reconstitution Method:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Stability:	12 months from date of despatch
Synonyms:	C7orf1; CPE-R2; CPETR2; HRVP1; RVP1
Summary:	Tight junctions represent one mode of cell-to-cell adhesion in epithelial or endothelial cell sheets, forming continuous seals around cells and serving as a physical barrier to prevent solutes and water from passing freely through the paracellular space. These junctions are comprised of sets of continuous networking strands in the outwardly facing cytoplasmic leaflet, with complementary grooves in the inwardly facing extracytoplasmic leaflet. The protein encoded by this intronless gene, a member of the claudin family, is an integral membrane protein and a component of tight junction strands. It is also a low-affinity receptor for Clostridium perfringens enterotoxin, and shares aa sequence similarity with a putative apoptosis-related protein found in rat. [provided by RefSeq, Jul 2008]



[View online »](#)