

Product datasheet for **TP724540**

Mouse TSLP Protein, hFc Tag

Product data:

Product Type:	Recombinant Proteins
Description:	Mouse TSLP Protein, hFc Tag
Expression Host:	HEK293
Tag:	C-Human Fc
Predicted MW:	The protein has a predicted molecular mass of 40.1 kDa after removal of the signal peptide. The apparent molecular mass of mTSLP-hFc is approximately 35-55 kDa due to glycosylation.
Purity:	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Reconstitution Method:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Stability:	12 months from date of despatch
Summary:	Enables cytokine activity. Predicted to be involved in several processes, including defense response to other organism; positive regulation of cytokine production; and positive regulation of signal transduction. Predicted to be located in extracellular region. Is expressed in ductus deferens; epididymis; ileum; and testis. Used to study atopic dermatitis; autoimmune hemolytic anemia; cryoglobulinemia; and membranoproliferative glomerulonephritis. Orthologous to human TSLP (thymic stromal lymphopoietin). [provided by Alliance of Genome Resources, Apr 2022]



[View online »](#)