

Product datasheet for **TP724630**

IL2 Receptor beta (IL2RB) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Interleukin-2 Receptor Subunit Beta/IL2RB (N-Avi&His)
Species:	Human
Expression cDNA Clone or AA Sequence:	A27-T240
Tag:	N-Avi&6His
Buffer:	Lyophilized from a 0.2 um filtered solution of PBS, pH7.4
Note:	Recombinant human interleukin-2 receptor subunit beta is produced by Human 293 Cells. The target gene encoding A27-T240 is expressed with a Avi&6His tag at the N terminus.
Stability:	12 months from date of despatch
Locus ID:	3560
UniProt ID:	P14784
Summary:	Interleukin-2 receptor subunit beta (IL2RB) belongs to the type I cytokine receptor family and it is expressed in mononuclear cells. IL2RB contains a fibronectin type-III domain. IL2RB / CD122 is involved in T cell-mediated immune responses and is present in 3 forms with respect to ability to bind interleukin 2. The low affinity form is a monomer of the alpha subunit and is not involved in signal transduction. The intermediate affinity form consists of an alpha/beta subunit heterodimer, while the high affinity form consists of an alpha/beta/gamma subunit heterotrimer. IL2RB is the receptor for interleukin-2. This beta subunit is involved in receptor mediated endocytosis and transduces the mitogenic signals of IL2. IL2RB might also be in association with IL15RA, involved in the stimulation of neutrophil phagocytosis by IL15.



[View online »](#)