

Product datasheet for **TP724714**

Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human SIGLEC-5 (C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	Glu17-Thr434
Tag:	C-His
Buffer:	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Note:	Recombinant Human SIGLEC-5 (C-His) is produced by Human Cells. The target gene encoding Glu17-Thr434 is expressed with a C-His tag.
Stability:	12 months from date of despatch
Synonyms:	SIGLEC5,OB-BP2,CD170,CD33L2



[View online »](#)