

Product datasheet for **TP724723**

Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human SIGLEC-10 (C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	Met 17-Thr 546
Tag:	C-His
Buffer:	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Note:	Recombinant Human SIGLEC-10 (C-His) is produced by Human Cells. The target gene encoding Met 17-Thr 546 is expressed with a C-His tag.
Storage:	Lyophilized protein should be stored at < -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Stability:	12 months from date of despatch
Synonyms:	SIGLEC10,MGC126774,PRO940,Siglec10,SLG2
Summary:	The siglecs (sialic acid-binding Ig-like lectins) are a distinct subset of the Ig superfamily with adhesion-molecule-like structure. We describe here a novel member of the siglec protein family that shares a similar structure including five Ig-like domains, a transmembrane domain, and a cytoplasmic tail containing two ITIM-signaling motifs. Siglec-10 was identified through database mining of an asthmatic eosinophil EST library. The Siglec-10-VAP-1 interaction seems to mediate lymphocyte adhesion to endothelium and has the potential to modify the inflammatory microenvironment via the enzymatic end products.



[View online »](#)