

Product datasheet for **TP724725**

Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Mouse SIGLEC-10 (C-mFc)
Species:	Mouse
Expression cDNA Clone or AA Sequence:	Gly17-Val546
Tag:	C-Mouse Fc
Buffer:	Lyophilized from a 0.2 μ m filtered solution of PBS, pH 7.4.
Note:	Recombinant Mouse SIGLEC-10 (C-mFc) is produced by Human Cells. The target gene encoding Gly17-Val546 is expressed with a C-mFc tag.
Stability:	12 months from date of despatch
Summary:	The siglecs (sialic acid-binding Ig-like lectins) are a distinct subset of the Ig superfamily with adhesion-molecule-like structure. We describe here a novel member of the siglec protein family that shares a similar structure including five Ig-like domains, a transmembrane domain, and a cytoplasmic tail containing two ITIM-signaling motifs. Siglec-10 was identified through database mining of an asthmatic eosinophil EST library. The Siglec-10-VAP-1 interaction seems to mediate lymphocyte adhesion to endothelium and has the potential to modify the inflammatory microenvironment via the enzymatic end products.



[View online »](#)