

## Product datasheet for **TP724738**

### Cynomolgus Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Cynomolgus MELTF (N-hFc)
Species:	Cynomolgus
Expression cDNA Clone or AA Sequence:	Gly111-Ala803
Tag:	N-Human Fc
Buffer:	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Note:	Recombinant Cynomolgus MELTF (N-hFc) is produced by Human Cells. The target gene encoding Gly111-Ala803 is expressed with a N-hFc tag.
Stability:	12 months from date of despatch



[View online »](#)