

Product datasheet for **TP724840**

Phospholipase A2 (PLA2G4A) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinat Human Cytosolic phospholipase A2(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-V138
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinat Human Cytosolic phospholipase A2 is produced by E.coli expression system and the target gene encoding S2-V138 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	5321
UniProt ID:	<u>P47712</u>



[View online »](#)