

## Product datasheet for **TP724841**

### FKBP51 (FKBP5) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human FKBP51 Protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	T2-E313
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human FKBP51 Protein is produced by E.coli expression system and the target gene encoding T2-E313 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	2289
UniProt ID:	<a href="#">Q13451</a>
Synonyms:	Peptidyl-prolyl cis-trans isomerase FKBP5; PPIase FKBP5; 51 kDa FKBP; FKBP-51; 54 kDa progesterone receptor-associated immunophilin; Androgen-regulated protein 6; FF1 antigen; FK506-binding protein 5
Protein Families:	Druggable Genome



[View online »](#)