

Product datasheet for **TP724844**

WASP (WAS) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Actin nucleation-promoting factor WAS(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	D121-E311
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Actin nucleation-promoting factor WAS is produced by E.coli expression system and the target gene encoding D121-E311 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	7454
UniProt ID:	P42768



[View online »](#)