

Product datasheet for **TP724847**

HSD11B1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human 11 β -HSD1 (C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	E25-A250
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human 11 β -HSD1 is produced by E.coli expression system and the target gene encoding E25-A250 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	3290
UniProt ID:	P28845



[View online »](#)