

## Product datasheet for **TP724857**

### EIF4G1 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human EIF4G1 Protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	L1401-N1599
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human EIF4G1 Protein is produced by E.coli expression system and the target gene encoding L1401-N1599 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	1981
UniProt ID:	<a href="#">Q04637</a>
Synonyms:	Eukaryotic translation initiation factor 4 gamma 1; eIF-4-gamma 1; eIF-4G 1; eIF-4G1; p220; EIF4F; EIF4G; EIF4GI
Protein Pathways:	Viral myocarditis



[View online »](#)