

## Product datasheet for **TP724861**

### TXNDC5 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human TXNDC5 Protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	E110-L432
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human TXNDC5 Protein is produced by E.coli expression system and the target gene encoding E110-L432 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	81567
UniProt ID:	<a href="#">Q8NBS9</a>
Synonyms:	Thioredoxin domain-containing protein 5; Endoplasmic reticulum resident protein 46; Thioredoxin-like protein p46; TLP46
Protein Families:	Druggable Genome



[View online »](#)