

## Product datasheet for **TP724862**

### Syntaxin 8 (STX8) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Syntaxin 8(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-G215
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Syntaxin 8 is produced by E.coli expression system and the target gene encoding A2-G215 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	9482
UniProt ID:	<a href="#">Q9UNK0</a>
Synonyms:	STX8
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	SNARE interactions in vesicular transport



[View online »](#)