

Product datasheet for **TP724864**

WASF3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Actin-binding protein WASF3(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	P2-K180
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Actin-binding protein WASF3 is produced by E.coli expression system and the target gene encoding P2-K180 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10810
UniProt ID:	Q9UPY6
Synonyms:	Protein WAVE-3; Verprolin homology domain-containing protein 3; Wiskott-Aldrich syndrome protein family member 3; WASP family protein member 3; WAVE3
Protein Families:	Druggable Genome
Protein Pathways:	Adherens junction, Fc gamma R-mediated phagocytosis



[View online »](#)