

Product datasheet for **TP724866**

NYREN18 (NUB1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human NEDD8 ultimate buster 1/NUB1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	R513-N615
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human NEDD8 ultimate buster 1/NUB1 is produced by E.coli expression system and the target gene encoding R513-N615 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	51667
UniProt ID:	Q9Y5A7
Synonyms:	Negative regulator of ubiquitin-like proteins 1; Renal carcinoma antigen NY-REN-18; NYREN18
Protein Families:	Druggable Genome



[View online »](#)