

Product datasheet for **TP724867**

LARG (ARHGEF12) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human ARHGEF12/LARG(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	I154-L575
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human ARHGEF12/LARG is produced by E.coli expression system and the target gene encoding I154-L575 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	23365
UniProt ID:	Q9NZN5
Synonyms:	Rho guanine nucleotide exchange factor 12; Leukemia-associated RhoGEF
Protein Pathways:	Axon guidance, Regulation of actin cytoskeleton, Vascular smooth muscle contraction



[View online »](#)