

## Product datasheet for **TP724873**

### Jaw1 (LRMP) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human IRAG2/LRMP Protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	Q143-P292
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human IRAG2/LRMP Protein is produced by E.coli expression system and the target gene encoding Q143-P292 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	4033
UniProt ID:	<a href="#">Q12912</a>
Synonyms:	Inositol 1,4,5-triphosphate receptor associated 2; Lymphoid-restricted membrane protein; Protein Jaw1; JAW1
Protein Families:	Transmembrane



[View online »](#)