

Product datasheet for **TP724875**

MAP4K5 Human Recombinant Protein

Product data:

| | |
|---------------------------------------|--|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant Human MAP4K5 Protein(C-His) |
| Species: | Human |
| Expression cDNA Clone or AA Sequence: | K395-N484 |
| Tag: | C-His |
| Buffer: | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4 |
| Note: | MADE TO ORDER Recombinant Human MAP4K5 Protein is produced by E.coli expression system and the target gene encoding K395-N484 is expressed with a His tag at the C-terminus. |
| Storage: | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles. |
| Stability: | 12 months from date of despatch |
| Locus ID: | 11183 |
| UniProt ID: | Q9Y4K4 |
| Synonyms: | Mitogen-activated protein kinase kinase kinase kinase 5; Kinase homologous to SPS1/STE20; MAPK/ERK kinase kinase kinase 5; KHS; MEK kinase kinase 5; MEKKK 5 |
| Protein Families: | Druggable Genome, Protein Kinase |



[View online »](#)