

Product datasheet for **TP724878**

Miz1 (ZBTB17) Human Recombinant Protein

Product data:

| | |
|---------------------------------------|---|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant Human Miz-1 Protein(C-His) |
| Species: | Human |
| Expression cDNA Clone or AA Sequence: | D2-A188 |
| Tag: | C-His |
| Buffer: | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4 |
| Note: | MADE TO ORDER Recombinant Human Miz-1 Protein is produced by E.coli expression system and the target gene encoding D2-A188 is expressed with a His tag at the C-terminus. |
| Stability: | 12 months from date of despatch |
| Locus ID: | 7709 |
| UniProt ID: | Q13105 |



[View online »](#)