

Product datasheet for **TP724881**

TRIM5 alpha (TRIM5) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human TRIM5 Protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-Q248
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombiant Human TRIM5 Protein is produced by E.coli expression system and the target gene encoding A2-Q248 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	85363
UniProt ID:	Q9C035
Synonyms:	Tripartite motif-containing protein 5; RING finger protein 88; RING-type E3 ubiquitin transferase TRIM5; RNF88
Protein Families:	Druggable Genome



[View online »](#)