

Product datasheet for **TP724882**

EIF4B Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human eIF4B Protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-R200
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombiant Human eIF4B Protein is produced by E.coli expression system and the target gene encoding A2-R200 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	1975
UniProt ID:	P23588
Synonyms:	Eukaryotic translation initiation factor 4B; eIF-4B
Protein Pathways:	mTOR signaling pathway



[View online »](#)