

## Product datasheet for **TP724885**

### LY75 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Lymphocyte antigen 75(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	C216-C501
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Lymphocyte antigen 75 is produced by E.coli expression system and the target gene encoding C216-C501 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	4065
UniProt ID:	<a href="#">O60449</a>
Synonyms:	C-type lectin domain family 13 member B; DEC-205; gp200-MR6; CD205; LY75; CD205, CLEC13B
Protein Families:	Transmembrane



[View online »](#)