

## Product datasheet for **TP724889**

### SAP18 Human Recombinant Protein

#### Product data:

|                                       |  |
|---------------------------------------|--|
| Product Type:                         | Recombinant Proteins   |
| Description:                          | Recombinant Human SAP18 Protein  |
| Species:                              | Human  |
| Expression cDNA Clone or AA Sequence: | M1-Y153  |
| Buffer:                               | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4  |
| Note:                                 | MADE TO ORDER   Recombinant Human SAP18 Protein is produced by E.coli expression system and the target gene encoding M1-Y153 is expressed.   |
| Storage:                              | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles. |
| Stability:                            | 12 months from date of despatch  |
| Locus ID:                             | 10284  |
| UniProt ID:                           | <a href="#">O00422</a>   |
| Synonyms:                             | Histone deacetylase complex subunit SAP18; 18 kDa Sin3-associated polypeptide; 2HOR0202; Cell growth-inhibiting gene 38 protein; Sin3-associated polypeptide p18; GIG38                                      |
| Protein Families:                     | Druggable Genome, Transcription Factors  |



[View online »](#)