

Product datasheet for **TP724890**

Pirin (PIR) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Pirin
Species:	Human
Expression cDNA Clone or AA Sequence:	K46-N230
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Pirin is produced by E.coli expression system and the target gene encoding K46-N230 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	8544
UniProt ID:	O00625
Synonyms:	Probable quercetin 2,3-dioxygenase PIR; Probable quercetinase; PIR
Protein Families:	Transcription Factors



[View online »](#)