

Product datasheet for **TP724893**

KIF1C Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Kinesin-like protein KIF1C
Species:	Human
Expression cDNA Clone or AA Sequence:	452-718
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Kinesin-like protein KIF1C is produced by E.coli expression system and the target gene encoding 452-718 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10749
UniProt ID:	O43896
Synonyms:	KIF1C
Protein Families:	Druggable Genome



[View online »](#)