

## Product datasheet for **TP724894**

### CBR3 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Carbonyl reductase [NADPH] 3
Species:	Human
Expression cDNA Clone or AA Sequence:	C150-T250
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Carbonyl reductase [NADPH] 3 is produced by E.coli expression system and the target gene encoding C150-T250 is expressed.
Stability:	12 months from date of despatch
Locus ID:	874
UniProt ID:	<a href="#">O75828</a>



[View online »](#)