

Product datasheet for **TP724895**

TRIO Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Triple functional domain protein
Species:	Human
Expression cDNA Clone or AA Sequence:	M2556-N2678
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Triple functional domain protein is produced by E.coli expression system and the target gene encoding M2556-N2678 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	7204
UniProt ID:	O75962
Synonyms:	PTPRF-interacting protein; TRIO
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)