

## Product datasheet for **TP724897**

### UBC Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Polyubiquitin-C
Species:	Human
Expression cDNA Clone or AA Sequence:	K559-V685
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Polyubiquitin-C is produced by E.coli expression system and the target gene encoding K559-V685 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	7316
UniProt ID:	<a href="#">P0CG48</a>
Synonyms:	UBC; Ubiquitin
Protein Families:	Druggable Genome
Protein Pathways:	PPAR signaling pathway



[View online »](#)