

Product datasheet for **TP724899**

CSPS (SULT1A3) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Sulfotransferase 1A3
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-L295
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Sulfotransferase 1A3 is produced by E.coli expression system and the target gene encoding M1-L295 is expressed.
Stability:	12 months from date of despatch
Locus ID:	6818
UniProt ID:	P50224



[View online »](#)