

## Product datasheet for **TP724900**

### CSPS (SULT1A3) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Sulfotransferase 1A4
Species:	Human
Expression cDNA Clone or AA Sequence:	1-295
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Sulfotransferase 1A4 is produced by E.coli expression system and the target gene encoding 1-295 is expressed.
Stability:	12 months from date of despatch
Locus ID:	6818
UniProt ID:	<a href="#">P50224</a>



[View online »](#)