

Product datasheet for **TP724901**

SAT1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Diamine acetyltransferase 1
Species:	Human
Expression cDNA Clone or AA Sequence:	1-171
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Diamine acetyltransferase 1 is produced by E.coli expression system and the target gene encoding 1-171 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	6303
UniProt ID:	<u>P21673</u>



[View online »](#)