

Product datasheet for **TP724905**

Ribosomal protein S10 (RPS10) Human Recombinant Protein

Product data:

| | |
|---------------------------------------|--|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant Human 40S ribosomal protein S10/ RPS10 |
| Species: | Human |
| Expression cDNA Clone or AA Sequence: | M1-Q165 |
| Buffer: | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4 |
| Note: | MADE TO ORDER Recombinant Human 40S ribosomal protein S10/ RPS10 is produced by E.coli expression system and the target gene encoding M1-Q165 is expressed. |
| Storage: | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles. |
| Stability: | 12 months from date of despatch |
| Locus ID: | 6204 |
| UniProt ID: | P46783 |
| Synonyms: | Small ribosomal subunit protein eS10 |
| Protein Pathways: | Ribosome |



[View online »](#)