

Product datasheet for **TP724906**

HAAO Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human 3-hydroxyanthranilate 3,4-dioxygenase/HAAO
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-G286
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human 3-hydroxyanthranilate 3,4-dioxygenase/HAAO is produced by E.coli expression system and the target gene encoding M1-G286 is expressed.
Stability:	12 months from date of despatch
Locus ID:	23498
UniProt ID:	P46952



[View online »](#)