

Product datasheet for **TP724907**

FXR1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human RNA-binding protein FXR1
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-Y353
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human RNA-binding protein FXR1 is produced by E.coli expression system and the target gene encoding M1-Y353 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	8087
UniProt ID:	P51114
Synonyms:	FMR1 autosomal homolog 1; hFXR1p; Fragile X mental retardation syndrome-related protein 1
Protein Families:	Druggable Genome



[View online »](#)