

## Product datasheet for **TP724913**

### BPTF Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human BPTF Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	E459-P657
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human BPTF Protein is produced by E.coli expression system and the target gene encoding E459-P657 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	2186
UniProt ID:	<a href="#">Q12830</a>
Synonyms:	Nucleosome-remodeling factor subunit BPTF; Bromodomain and PHD finger-containing transcription factor; Fetal Alz-50 clone 1 protein; Fetal Alzheimer antigen; FAC1; FALZ
Protein Families:	Druggable Genome



[View online »](#)