

Product datasheet for **TP724916**

NFATC4 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human NF-ATc4/NFATc4
Species:	Human
Expression cDNA Clone or AA Sequence:	G101-W225
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human NF-ATc4/NFATc4 is produced by E.coli expression system and the target gene encoding G101-W225 is expressed.
Stability:	12 months from date of despatch
Locus ID:	4776
UniProt ID:	Q14934



[View online »](#)