

## Product datasheet for **TP724919**

### MAPKAP Kinase 3 (MAPKAPK3) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human MAP kinase-activated protein kinase 3
Species:	Human
Expression cDNA Clone or AA Sequence:	K179-Q382
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human MAP kinase-activated protein kinase 3 is produced by E.coli expression system and the target gene encoding K179-Q382 is expressed.
Stability:	12 months from date of despatch
Locus ID:	7867
UniProt ID:	<a href="#">Q16644</a>



[View online »](#)