

## Product datasheet for **TP724921**

### **HABP4 Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human Intracellular hyaluronan-binding protein 4
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	M1-K171
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	MADE TO ORDER   Recombinant Human Intracellular hyaluronan-binding protein 4 is produced by E.coli expression system and the target gene encoding M1-K171 is expressed.
<b>Storage:</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	22927
<b>UniProt ID:</b>	<a href="#">Q5JVS0</a>
<b>Synonyms:</b>	IHABP-4; IHABP4; Hyaluronan-binding protein 4 ; Ki-1/57 intracellular antigen; HABP4



[View online »](#)